



# Repeatability of Coulter Z2 and Vi-CELL Counts Using ViaCheck™ 0.5e+6 Concentration Controls

## INTRODUCTION


ViaCheck™ Viability and Concentration Controls (Bangs Laboratories) are available in a range of viability and concentration levels, the latter of which include 1e+6 beads/mL, 4e+6 beads/mL and 8e+6 beads/mL standards. Before introducing a new, lower-count standard (0.5e+6 beads / mL) to the ViaCheck product line, we endeavored to both confirm the repeatability of lower-end counting capabilities of the instruments featured in our in-house Quality Assurance program, and ensure that instrument performance supported intended product specifications.

## METHOD AND RESULTS

All samples were prepared and run by operators trained in proper sample preparation (including mixing and pipetting), instrument set-up and instrument operation<sup>3,4</sup>. Particle / cell counting instruments, including the Z2 (Coulter) and Vi-CELL® XR (Coulter) were properly calibrated and configured using standardized instrument settings<sup>1,6</sup>.

## REPEATABILITY: INTRA-INSTRUMENT PRECISION


To determine the range of within-instrument variation for the new lower-count ViaCheck standard (0.5e+6 beads/mL (VC50N), a total of 102 samples were run (51 on the Vi-CELL XR and 51 on the Z2). Precision was found to be acceptable for both the Z2 (<2%)<sup>2</sup> and Vi-CELL XR (<6%)<sup>5</sup>.




RESULTS	Z2	Vi-CELL® XR
Avg conc. (n = 51 [each])	5.32e+05 (beads/mL)	5.36e+05 (beads/mL)
Std. Dev.	9.90e+03 (beads/mL)	2.21e+04 (beads/mL)
%CV	1.86%	4.11%

## CONFIRMATION OF MANUFACTURING PRODUCT SPECIFICATIONS

Runs were performed against a product (count) specification range of 0.41e+6 – 0.59e+6 beads/mL (Coulter Z2). Per in-house policy, QC runs for each product Lot include both formal counts on the Z2, and confirmatory checks on the Vi-CELL XR. Each Lot is dispensed under a Squeglia Zero-Based Acceptance Sampling Plan; dispensing runs for 20mL and <1mL SingleShot™ volumes of the new 0.5e+6 beads/mL ViaCheck Concentration Control were evaluated. Product (count) specifications were found to be acceptable for both product formats, 20mL bottles and SingleShot vials.



Manufacturing Specification	0.41e+6 – 0.59e+6 beads / mL	
RESULTS	Z2	Vi-CELL XR
Avg conc. (n = 34 [each])	5.20e+05 (beads/mL)	5.28e+05 (beads/mL)
Std. Dev.	5.79e+03 (beads/mL)	2.20e+04 (beads/mL)
%CV	1.11%	4.17%



Manufacturing Specification	0.41e+6 – 0.59e+6 beads / mL	
RESULTS	Z2	Vi-CELL XR
Avg conc. (n = 13 [each])	5.10e+05 (beads/mL)	5.34e+05 (beads/mL)
Std. Dev.	5.97e+03 (beads/mL)	2.47e+04 (beads/mL)
%CV	1.17% (beads/mL)	4.62% (beads/mL)

## CONCLUSIONS

Both the Z2 and Vi-CELL XR demonstrated acceptable precision for counts of 0.5e+6 beads/mL, with the Z2 exhibiting slightly less variation. This performance difference is as documented in the literature, and is likely due to the underlying instrument methodology (Z2 – impedance; Vi-CELL XR – imaging). The additional reason for the difference in variation is likely and the number of counts typically performed in a run (Z2 – rate measured in kilocounts/sec, <100,000 counts; Vi-CELL – 50 fields/images, <2000 counts). All results fell within the specified count range for manufacturing (0.41e+6 – 0.59e+6 beads/mL), confirming the acceptability of the specification.

## REFERENCES & FURTHER READING

1. Beckman Coulter, Inc. (2011) *Vi-CELL XR Cell Viability Analyzer Reference Manual* (Publication No. 383674BA). Brea, CA
2. Beckman Coulter, Inc. (2007) *PROService™ Advisor: Z-Series* (Publication No. 05-2007-045). Brea, CA.
3. Bangs Laboratories, Inc. (2017) *ViaCheck for Cell Viability Analyzers: Best Practices*. (Tech Support Doc 0711) Fishers, IN.
4. Bangs Laboratories, Inc. (2017) *Handling & Pipetting Concentration Standards*. (Tech Support Doc 0706) Fishers, IN.
5. Lew C, Gomez JA, Rhyner MN. (2012) *Instrument-to-instrument variability in the Vi-CELL automated viability analyzer*. (Publication No. IB-17279A) Brea, CA. Beckman Coulter, Inc.
6. Kilbride K, Anglea B, Bavender A. (2017) *Optimization of Vi-CELL® XR settings for calibration checks using ViaCheck™ Controls*. (Application Note 0708.) Fishers, IN. Bangs Laboratories, Inc.